USER MANUAL

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INSTRUCTION MANUAL

MODE D'EMPLOI

MANUAL DE USARIO

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★ TriggerFairy®★

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★ TriggerFairy®★

USER MANUAL

2.1 Upper neck – rear of neck	Page 17
2.2 Lower neck to scapular spine	Page 19
2.3 Side of the jaw	Page 21
2.4 Shoulder blade	Page 22
2.5 Between the shoulder blades	Page 23
2.6 Below the shoulder blade	Page 25
2.7 Outer side of chest	Page 26
2.8 Middle of back	Page 27
2.9 Under the collarbone	Page 28
General remarks	Page 29

2.1. Massage zone on upper neck – rear of neck

The following massage techniques can be used:

- ★ Ischemic compression
- ★ Precise massage strokes
- ★ Pressure and movement technique

Area massaged:

You will be working on the muscles of your entire cervical spine and upper thoracic spine.


Feel your way along this area of your spine so that you get to know it better. Place one or two fingers on the back of your head and nod your head a few times.

Every time you nod your head downwards, you'll be able to feel a cord-like structure protruding. (1)

Place your fingers directly beside this 'cord' and move slowly down your neck.

Feel the muscles located there by applying pressure from the side as well as from behind. (2,3)

You will be massaging this entire area down to your lower neck area, where you can feel the mostly very firm and thick part of the trapezius.

This is the area that most people instinctively knead and massage when their neck hurts. (4)

Caution:

Start the massage gently, especially if you have never done self-massage before or suffer from sensitive muscles. Avoid pressing directly on the spine.

Massage with the TriggerFairy®


Hold the Trigger Fairy® in both hands.

Place one hand on the bow. Place the index finger of this hand on the outside or top of the Trigger Fairy®. Your other hand should be on the handle. (5)

Hold the Trigger Fairy® lightly without gripping too tightly, which would lead to tension in your hands and arms.

Perform the following techniques applying pressure from the side as well as pressure from behind. (6, 7, 8, 9)

Ischemic compression

Apply pressure with the head of the TriggerFairy® to the desired area and maintain the pressure for 30–60 seconds. Breathe calmly and try to relax. You can also make small, slow circular movements by moving the TriggerFairy® handle.

Precise massage strokes

Use the TriggerFairy® to make slow strokes from top to bottom of the massage area. Use the index finger of the massaging hand to exert pressure on the downward stroke.

Slow down your stroke if you come to a painful spot.

Pressure and movement technique

Exert pressure in the desired area and perform the following movements.

Turn your head slowly to the left and right. (10)

Nod your head slowly backwards and forwards. (11)

Experiment with the range of your head movements. (12)

2.2 Massage zone from lower neck to scapular spine

The following massage techniques can be used:

- ★ Pressure and movement technique
- ★ Ischemic compression
- ★ Precise rolling movements

Area massaged:

The massage area is on the shoulder and lower neck.

You will be massaging the area next to the superior angle of scapula (angulus superior), or medial angle, and above the scapular spine.

To feel the scapular spine, place your hand over the opposite shoulder and use your fingers to feel the bony line running more or less horizontally on your shoulder blade. (13)

To feel the superior angle of scapula, move along the bone towards the middle of the body and then move upwards a little.

There you can feel a bony node-like structure. This is the superior angle of your shoulder blade. (14)

In feeling for these anatomical landmarks, proceed slowly and exert only gentle pressure. Too much pressure and haste will prevent you from feeling them properly.


General massage tips:

This area is often very sensitive, especially if you are sitting at a desk for long periods or experience a lot of stress in daily activity. Therefore make sure you breathe gently and evenly during the massage and keep your shoulders loose.

Massage with the TriggerFairy®

Hold the TriggerFairy® with your opposite hand and make sure that the arm on the side to be massaged is hanging loose and free.

Place the TriggerFairy® on the desired area and perform one of the following techniques. (15)

Ischemic compression

Press the head of the Fee onto the desired area and exert pressure. Maintain pressure for 30–60 seconds.

Breathe deeply and evenly. With each breath, try to relax your shoulder/neck a little more and allow them to sink.

Pressure and movement technique

This technique can be performed with the TriggerFairy® in two different positions.

Do it once using pressure from above and once using pressure from the front.

Obviously other muscle fibres will also be massaged in the process. To make sure you have attended to all the muscles in this area, you should carry out both versions at least once.

This technique is usually accompanied by shoulder movements. To this end you can ...

- ★ Rotate the shoulders.
- ★ Lift and lower your shoulders.
- ★ Raise your arm out to the side and lower it again.


Carry out all these movements very slowly. Make them small, large and every stage inbetween.

Pressure from above

Hold the TriggerFairy® with your opposite hand and ensure the arm on the side to be massaged is hanging loose and free.

In the area of your superior angular of scapula, rotate, lift and lower your shoulder. (15)

In the area directly over the scapular spine, lift and lower your arm.

Pressure from the front

Hold the TriggerFairy® with your opposite hand. Hold it loosely on the bow and allow the handle to rest on your stomach or chest. (16)

Exert pressure on the rod in the area desired and rotate your shoulder. (17)

Precise massage strokes

With this technique you can cover the entire upper shoulder area. Perform the massage on bare skin so that the rod can glide smoothly.

Using your opposite hand, hold the TriggerFairy® near to your head. Place your index finger on the upper side of the TriggerFairy® and place your thumb on the outside.

Place the TriggerFairy® over your superior angle of scapula. (18)

Make long, slow massage strokes along your neck. Vary the direction of pressure, applying it from the front, from above and from behind.

2.3 Massage zone on the side of the jaw

The following massage techniques can be used:

- ★ Ischemic compression
- ★ Precise massage strokes
- ★ Pressure movement technique

Area massaged:

You will be working on the area on the side of your jaw and the chewing muscles located there. To get to know this area and its muscles better, feel it while it is tensed.

Place your finger lightly on your jaw and clamp your teeth together several times.

With every bite, you'll be able to feel your jaw muscles tense under your fingers. Try to feel the movement and boundaries of these muscles. In the techniques that follow, you will be exploring and massaging this entire area.

Massage with the TriggerFairy®

Hold the Fee in both hands, with one hand on the bow and the other on the handle.

The hand on the bow should be gripping with the thumb underneath and the index finger on top. (19)

19

Pressure and movement technique

Apply increasing pressure in the desired area and perform the following movements.

Open and close your mouth and clamp your teeth together. (20)

Open your jaw and move it left and right. (21)

20

21

Place the head of the Fee on the desired area and perform one of the following techniques.

Ischemic compression

Using the head of the Fee, exert pressure on the desired area and maintain it for 30–60 seconds.

Breathe calmly and try to relax.

You can also carry out small, slow circular movements by moving the handle of the Fee.

Precise massage strokes

Move the TriggerFairy® slowly over the massage area from top to bottom.

Use the index finger of the massaging hand to exert pressure on the downward stroke. Slow down your strokes when you come across a sensitive area.

2.4 Massage zone Shoulder blade

Area massaged:

Place the fingertips of your opposite hand over your shoulder so that you can feel the scapular spine. This is the bony ridge running more or less horizontally along your shoulder blade. (13)

Massage the area below this ridge.

Massage with the TriggerFairy®

Hold the TriggerFairy® with your opposite hand. (22)


Place the TriggerFairy® in the desired area and perform one of the following techniques.

Ischemic compression

Using the head of the TriggerFairy®, exert pressure on the desired area and maintain it for 30–60 seconds. (23)

Breathe calmly and try to relax. You can also carry out small, slow circular movements by moving the handle of the Fee.


Pressure and movement technique

Apply increasing pressure to the desired area and perform the following movements.

Rotate your upper arm as far inwards and as far outwards as possible. (24+25)

Raise and lower your arm. (26)

Note:

Many people rotate only their lower arm. This has no effect on the shoulder and therefore no effect on your shoulder blade's musculature.

If you ensure you rotate the elbow, you will also be rotating your upper arm. Alternatively you can raise and lower your arm.

2.5 Massage zone between the shoulder blades

The following massage techniques can be used:

- ★ Pressure and movement technique
- ★ Ischemic compression

Area massaged:

This massage area is located between the spine and the shoulder blade of the side in question. Work the area from the superior angle of scapula down to below the lower angle of scapula. (14+27)

Do not massage your spine directly. Start right next to it.

Massage with the TriggerFairy®

Hold the TriggerFairy® with your opposite hand.
Place the head of the TriggerFairy® on the desired area. (28)

Perform one of the following techniques.

Ischemic compression

Using the head of the TriggerFairy®, exert pressure on the desired area and maintain it for 30–60 seconds.

Breathe calmly and try to relax.

You can also carry out small, slow circular movements by moving the handle of the Fee.

Pressure and movement technique

Apply increasing pressure to the desired area and perform the following movements.

Perform circular movements with your shoulder. (30, 31, 32)

Move your shoulder blade by moving your arm backwards and forwards. (33)

2.6 Massage zone below the shoulder blade

The following massage techniques can be used:

- ★ Ischemic compression

Area massaged:

You will be massaging the area below your shoulderblade, particularly along its entire inner edge (margo medialis).

This extends from the superior to the lower angle of scapula.

To feel the superior angle of scapula (angulus superior), place one or two fingers of your opposite hand over your shoulder and try to feel the bony nodule-like structure located there. (14)

To feel the lower angle of scapula, reach behind your back with your opposite hand and try to feel the lowest point of your shoulder blade, its small 'corner'. (27)

If your shoulder is in the neutral position, you will not be able to reach this area below the shoulder blade. You therefore have to cross the arm on the side in question in front of your body and keep it there.

This arm movement forces the shoulder blade outwards and allows you to reach the area in question. (34)

General tips for massaging this area

A massage here can often be quite painful.

All the more important, therefore, to apply an appropriate amount of pressure when massaging, to breathe gently and evenly, and to relax your face and jaw muscles while performing the massage.

Massage with the TriggerFairy®

Hold the TriggerFairy® with your opposite hand. (33)

Place the head of the Trigger Fairy® in the desired area.

Perform one of the following techniques.

Ischemic compression

Use the head of the TriggerFairy® to apply pressure in the desired area and maintain the pressure for 30–60 seconds. (35)

You can also carry out small, slow circular movements by moving the handle of the TriggerFairy®.

Breathe calmly and try to relax.

With each out breath, try to relax your shoulder a little more.

2.7 Massage zone on outer side of chest

The following massage techniques can be used:

- ★ Ischemic compression
- ★ Pressure and movement technique

Area massaged:

As suggested in the heading to this section, you will be massaging the outer side of your chest.

More precisely, this is the area below as well as left and right of your coracoid process. This is the bony, nodule-like structure that you feel if you slowly trace a path with your fingers upwards from the lower outer side of the chest.

Massage the 10–15cm lying diagonally below.

General tips for massaging this area – and common mistakes

Allow the shoulder and arm of the side to be massaged to hang loose.

Breathe deeply and evenly into your chest during the massage.

Massage with the TriggerFairy®

Hold the TriggerFairy® with your opposite hand. Grip the TriggerFairy® at the bow and place your index finger on the bow. (36)

Perform one of the following techniques.

Ischemic compression

Use the head of the TriggerFairy® to apply pressure in the desired area and maintain the pressure for 30–60 seconds.


Breathe calmly and try to relax.

Pressure and movement technique

Apply increasing pressure to the desired area and perform the following movements.

Move your shoulder backwards several times. (37)

Raise your shoulders several times. (38)

2.8 Massage zone in middle of back

The following massage techniques can be used:

- ★ Ischemic compression
- ★ Pressure and movement technique

Area massaged:

You will be massaging your back muscles in the area from the lower angle of scapula into the lumbar vertebral area.

Your lower angle of scapula is the deepest point of your shoulder blade.

You can feel this bony point by reaching behind your back with your opposite hand to touch your shoulder blade and then moving your shoulder blade in small circles. (27)

Massage with the TriggerFairy®

Hold the TriggerFairy® with both hands.

Hold the handle with the hand opposite to the side to be massaged. (39)

Grip the bow with the other hand.

Place the TriggerFairy® in the desired area and perform one of the following techniques.

Ischemic compression

Use the head of the TriggerFairy® to apply pressure in the desired area and maintain the pressure for 30–60 seconds. (40)

Breathe calmly and try to relax.



You can also make small, slow circular movements by moving the handle of the TriggerFairy®.

Imagine that you are breathing directly into the area where the TriggerFairy® is situated.

Pressure and movement technique

Apply increasing pressure to the desired area and perform the following movements.

Bend backwards and forwards a few times. (41+42)

2.9 Massage zone under the collarbone

The following massage techniques can be used:

- ★ Precise massage strokes
- ★ Pressure and movement technique

Area massaged:

You will be working on the area directly below your collarbone.

To get a feel for this area, place your finger on the upper end of your sternum and feel the bone that runs from here to the side of the body. That is your collarbone. (43+44)


43

44

Now apply pressure with your fingertips to the area directly below the collarbone. (45)

45

Start next to the sternum. After a few centimetres you will come across a small bony nodule.

Massage the area directly below the collarbone and next to the sternum, going right up to the coracoid process.

Massage with the TriggerFairy®

Hold the Fee with your opposite hand roughly in the middle of the bow. Your index finger should be on the top of the TriggerFairy® and your thumb on the outer side. (46)

46

Place the head of the TriggerFairy® at a 45° angle to the massage area.

The handle of the rod should be resting on your upper body. Place the TriggerFairy® in the desired area and perform one of the following techniques.

Precise massage strokes

Use the rod to perform very small horizontal massage strokes of 1–2cm towards the middle of the body.

Make sure you simply ,shift' the skin rather than stroking it.

Pressure and movement technique

Perform small massage strokes (see above) and ...
... at the same time pull your shoulders slowly backwards.

Massage tips for this area

Make sure you keep the TriggerFairy® diagonal. This will help the massage greatly.

Don't tense your fingers too much when you grip the TriggerFairy®, as this can cause slight cramp in the hand. Rather press down on the Fee from above using a ,loose' hand.

Tips for using the TriggerFairy®

- ★ The TriggerFairy® is a wellness product. It should not be used as a substitute for professional or therapeutic treatment.
- ★ Use of the TriggerFairy® does not constitute any kind of curative treatment or therapy.
- ★ Do not exert pressure on your spine, bones or joints or other sensitive areas of the body.
- ★ Do not use while driving a vehicle or operating machinery.
- ★ Keep the TriggerFairy® out of the reach of children. Users must be aged 16 or over.
- ★ Do not use the TriggerFairy® if you have any illness affecting your muscles, joints or nerves or if you have any inflammation of the musculoskeletal system. If you suffer from any such condition but would still like to use the TriggerFairy®, consult your doctor in advance.
- ★ This device has been developed solely for the relaxation and exercise of the muscles.
- ★ Do not use for any unauthorised purpose.
- ★ Maximum tensile load 11 lbs.
- ★ Two years warranty.

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